



ANALIZA BIAŁEK



100006909

## Laboratory of Analysis

**Client:** Camo Chem  
**Order number:** 100006909  
**Name:** NOT MARKED  
**Number of samples:** 1 vial

**Batch number:** -----  
**Receive date:** 2025-01-28  
**Start date:** 2025-02-03  
**Enda date:** 2025-02-05

## Sample identification

Substance	Method of taking a sample	Form
GW-501516	Sample provided by customer	White powder in a glass vial


## Results

Substance	Sample name	Content	Purity
GW-501516	NOT MARKED	-----	98%

## Research comment

The test sample was applied in the amount of 2 $\mu$ l to the BioBasic 8 Thermo Scientific column. During the course of the analysis with the method a linear gradient in the solutions of H<sub>2</sub>O and Acetonitrile with the addition of 0.1% trifluoroacetic acid, a peak of the sought substance appeared. Its purity, when compared to the reference chromatogram, is 98%. Height the measurement error is + - 1%. The actual state of the sample sent before the test is presented on the photo below

## Photo of the test sample



**Password**

X83M29KV


**Analiza Białek sp. z o.o.**  
**KRS:** 0000966168  
**NIP:** 8992922970  
**REGON:** 521751875

**Headquarters and Laboratory:**  
ul. Duńska 9, 54-427 Wrocław  
Budynek DELTA

**Lab technician:**  
Analysis and research specialist  
mgr Arkadiusz Zajęc



## Chromatogram



## Method

Time	H2O + 0,1% TFA	ACN + 0,1% TFA
0min	60%	40%
2min	60%	40%
27min	5%	95%
30min	5%	95%

## Password

X83M29KV